

Silicone Adhesive Foam Dressings

Silicone Adhesive Foam Dressings: Highly conformable dressing for wounds with light to moderate exudate. Perforated, soft silicone gel adhesive provides gentle, but firm, adhesion and minimizes pain and trauma to the wound and skin upon removal. The dressing may be lifted and repositioned without losing its adherent properties. Breathable, yet moisture- and bacteria-resistant outer film layer seals the wound margins. The absorbent inner hydrocellular foam pad maintains a moist wound environment to help promote healing. McKesson silicone foam dressings are an excellent choice when fragile skin is a concern.



Indications: Wide range of acute and chronic wounds like venous leg ulcers, pressure ulcers, diabetic ulcers and traumatic wounds.

Contraindications: Patients with a known hypersensitivity to the product itself or its components.

Change Frequency: Change dressing as directed by a healthcare professional.

Item #3728987

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

McKesson ConnectTM
www.connect.mckesson.com
Customer Support
 855.625.7385 | Mon–Fri, 7am–9pm (CT)

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

Ordering Information

MFR #	US PHARMA ITEM #	DESCRIPTION	QTY	HCPCS	EXUDATE LEVEL*
4842	3728987	Hydrocellular Foam Dressings, Silicone Adhesive, With Border, Gentle Adhesive, Sterile, 3 in x 3 in	1/ea, 10/bx, 20bx/cs	A6212	●●○
4843	3729001	Hydrocellular Foam Dressings, Silicone Adhesive, With Border, Gentle Adhesive, Sterile, 4 in x 4 in	1/ea, 10/bx, 20bx/cs	A6212	●●○
4844	3729019	Hydrocellular Foam Dressings, Silicone Adhesive, With Border, Gentle Adhesive, Sterile, 6 in x 6 in	1/ea, 10/bx, 20bx/cs	A6213	●●○
4845	3729027	Hydrocellular Foam Dressings, Silicone Adhesive, With Border, Sacral, Gentle Adhesive, Sterile, 7 in x 7 in	1/ea, 10/bx, 20bx/cs	A6214	●●○
4862	3728771	Hydrocellular Foam Dressings, Silicone Adhesive, Without Border, Gentle Adhesive, Sterile, 3 in x 3 in	1/ea, 10/bx, 20bx/cs	A6209	●●○
4863	3728789	Hydrocellular Foam Dressings, Silicone Adhesive, Without Border, Gentle Adhesive, Sterile, 4 in x 4 in	1/ea, 10/bx, 20bx/cs	A6209	●●○
4864	3729068	Hydrocellular Foam Dressings, Silicone Adhesive, Without Border, Gentle Adhesive, Sterile, 6 in x 6 in	1/ea, 10/bx, 20bx/cs	A6210	●●○
4865	3728839	Hydrocellular Foam Dressings, Silicone Adhesive, Without Border, Sacral, Gentle Adhesive, Sterile, 7 in x 7 in	1/ea, 10/bx, 20bx/cs	A6211	●●○
4872	3728847	Lite Hydrocellular Foam Dressings, Silicone Adhesive, With Border, Gentle Adhesive, Sterile, 3 in x 3 in	1/ea, 10/bx, 20bx/cs	A6212	●○○
4873	3728870	Lite Hydrocellular Foam Dressings, Silicone Adhesive, With Border, Gentle Adhesive, Sterile, 4 in x 4 in	1/ea, 10/bx, 20bx/cs	A6212	●○○
4874	3728888	Lite Hydrocellular Foam Dressings, Silicone Adhesive, With Border, Gentle Adhesive, Sterile, 6 in x 6 in	1/ea, 10/bx, 20bx/cs	A6213	●○○
4892	3728896	Lite Hydrocellular Foam Dressings, Silicone Adhesive, Without Border, Gentle Adhesive, Sterile, 3 in x 3 in	1/ea, 10/bx, 20bx/cs	A6209	●○○
4893	3728904	Lite Hydrocellular Foam Dressings, Silicone Adhesive, Without Border, Gentle Adhesive, Sterile, 4 in x 4 in	1/ea, 10/bx, 20bx/cs	A6209	●○○
4894	3728920	Lite Hydrocellular Foam Dressings, Silicone Adhesive, Without Border, Gentle Adhesive, Sterile, 6 in x 6 in	1/ea, 10/bx, 20bx/cs	A6210	●○○

EXUDATE LEVEL KEY: ●○○ = Light ●●○ = Moderate ●●● = Heavy

*Exudate level indicates highest level of exudate that the dressing can accommodate

Silicone Adhesive Bandages

McKesson silicone adhesive bandages help protect cuts and scrapes by conforming to the skin with an absorbent, non-stick pad. The sterile, individually wrapped packages help promote cleanliness and further protect against infection.



Features and Benefits

- Features comfortable fit, long-lasting adhesive
- Flexible fabric bandages adhere securely, yet remove cleanly and are easy to reapply
- Transparent strips
- Not made with natural rubber latex

Item #3728979

Ordering Information

MFR #	US PHARMA ITEM #	DESCRIPTION	QTY
16-4840	3728979	Silicone Adhesive Bandages, Transparent, Strips, Sterile, 1 in x 3 in	1 ea, 10 ea/bx, 32 bx/cs

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

McKesson ConnectTM
www.connect.mckesson.com
Customer Support
 855.625.7385 | Mon–Fri, 7am–9pm (CT)


The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson or its affiliates (“McKesson”). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

Contact Layer Dressings

McKesson contact layer dressings are sterile, versatile dressings that can be used on a variety of different wounds when used with an outer absorbent dressing. The thin transparent dressing is made from polyurethane film and silicone gel with an open perforated structure. These dressings conform to the shape of the wound, while the open perforated structure allows exudate to pass vertically into the outer absorbent dressing and enables easy delivery of topical treatments. They prevent the outer absorbent dressing from sticking to the wound or the skin around the wound. They can also be used as a protective layer on superficial wounds that are not discharging and on areas of skin that are very fragile.



Features and Benefits

- Channels allow exudate to transfer vertically into a secondary absorbent dressing and aid delivery of topical treatments
- Designed to allow wound fluid to easily transfer into a secondary dressing
- Protects the wound and the skin
- Prevents an outer dressing from sticking to the wound
- Because the dressing is translucent, it allows the caregiver to assess the healing progress without having to remove the dressing
- Made without natural rubber latex

Item #3728987

Ordering Information

MFR #	US PHARMA ITEM #	DESCRIPTION	QTY
4816	3728953	Contact Layer Dressings One-Sided, Silicone Adhesive, Sterile, 3 in x 4 in	1 ea, 10 ea/bx, 4 bx/cs
4817	3728961	Contact Layer Dressings One-Sided, Silicone Adhesive, Sterile, 4 in x 7.2 in	1 ea, 10 ea/bx, 7 bx/cs

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

McKesson ConnectTM
www.connect.mckesson.com

Customer Support
 855.625.7385 | Mon–Fri, 7am–9pm (CT)

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.